

BATHTUB CRATE OPENING & UNPACKING INSTRUCTIONS

1. Open all 5 plywood plates by unfastening the fasteners with battery powered screwdriver. Start with plate marked with „**open this side first**”.
2. Release the decorative front panel from the adhesive tape. DO NOT USE knife or ANY other sharp instruments which could damage the front panel!
3. Carefully remove the decorative front panel from the crate
4. Use your hands to unfasten the 2 white bathtub feet which are tightened through the wooden planks (please refer to the picture below).
5. Tilt the crate so that the bathtub stands on its feet.
6. Loosen corner mounting plates located in all corners, by untightening the fasteners using battery powered screwdriver.
7. Carefully remove the bathtub from the crate
8. Screw on the 2 white bathtub feet which were previously removed according to p.4.
9. To install the bathtub, please, follow bathtub installation instructions.

